

Pipette-Adapter for Capillaries (PAC)

The **Pipette-Adapter for Capillaries (PAC)** assists in the targeted and dropwise draining of capillaries on POCT-test strips or into vessels.

It is suitable for plastic capillaries with various outer diameters, variable lengths as well as common preparations.

On request also available already completed with capillary.



| Item No. | Description | Note |
|----------|--|----------------|
| 061625 | Pipette-Adapter for Capillaries (PAC) Ø 1.3-1.4 mm | e.g. for EK 10 |
| 061630 | Pipette-Adapter for Capillaries (PAC) Ø 1.7-1.8 mm | e.g. for EK 20 |



Blood sampling and Pipetting:

(Insert the capillary up to the edge into the tip of the PAC before use.)

Prepare and disinfect the puncture site. Perform puncture with a safety lancet and discard the first drop of blood.

Hold the capillary horizontally for sampling so that it fills up by the capillary effect. Make sure that the mouth of the PAC is not covered during sampling!

After filling the capillary entirely, place your thumb gently on the upper mouth of the PAC so it is completely closed.

The (drop wise) draining is carried out by generating a slight gauge pressure created by the thumb.

Attention! To avoid the exposure of infectious blood with the operator, the PAC with the filled capillary should be held slightly downwards to avoid a withdrawal of the liquid inside the adapter.

While using the PAC the wearing of personal protective clothing (disposable gloves, etc.) is strongly recommended. For single-use only.

