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Instructions for use for KABEVETTE® Blood gas

apply to the products KABEVETTE® Blood gas BG 311 KA 2 individually sterile packed

Intended use

The product KABEVETTE® blood gas individually sterile packed is a sample container for receiving, transporting and, if necessary, storing a blood sample from the human body. It is a single-use IVD product, for the determination of blood gases and electrolytes, which is assigned to product class A according to (EU) 2017/746 Annex VIII Para. 2,5 (Rule 5c) and is placed on the market in a sterile condition.

The product is used in a professional environment by medically trained professionals. The use of the KABEVETTE® blood gas is intended for patients of all ages and genders, provided that the complete filling of the sample vessel with blood is ensured.

Features

The blood collection system KABEVETTE® Blood Gas BG 311 KA 2 individually sterile packed is a CE-marked, EO-sterilized, individually packed IVD product. The system consists of a tube with piston rod, each made of plastic, a Luer cone, a natural colored stopper with Luer cone and cone cap and is prepared with balanced heparin (12 - 30 IU/ml according to DIN EN ISO 6710:2017).

Notes

- Wear protective gloves and, if necessary, other protective clothing!
- Inspect blister packaging in advance and do not use tubes from damaged or previously accidentally opened blister packaging. These must be disposed of immediately.
- Use the product only within the indicated shelf life.
- Any violent handling of the product must be avoided.
- The product is intended for blood collection only.
- The used or contaminated cannula and accessories are to be disposed of in the appropriate disposal container for sharps.
- The product is for single use only! Multiple use of the product can lead to injuries, infections or death. The product must not be sterilized again.
- The exact filling quantity of the tube must be observed!

Application

Ensure that all necessary components of a blood collection are prepared.

The blood collection tube can be completely filled in different ways, e.g. by connecting the Luer cone of the KABEVETTE® Blood Gas BG 311 KA 2 with a Luer disposable cannula, with a Luer connection of the three-way stopcock of a permanent indwelling catheter already placed in the patient's vein, or using a custom-made Luer adapter with a cannula for multiple sampling.

When used in combination with a three-way stopcock, ensure that the tubing is completely filled with blood up to the connection point of the Luer cone of the KABEVETTE® Blood Gas. If the collection is performed using a wing cannula, the KABEVETTE® Blood Gas should not be the first blood collection tube to be taken due to the possible formation of air bubbles.

The following is an example of the application when using a disposable Luer cannula:

Prepare puncture site and do not touch after disinfection.

- 1) Pull the blister pack apart with your thumbs and forefingers at the point marked with an arrow and remove KABEVETTE®.
- 2) Remove the cone cap from the KABEVETTE® with a slight pulling and turning motion.
- 3) Attach a commercially available Luer cannula to the KABEVETTE® and remove the protective cannula cap.
- 4) Puncture the blood vessel with the ground edge of the cannula pointing upwards. During blood collection, hold the tube in one hand and gently pull the piston rod to the stop with the other hand.
- 5) Remove the tourniquet. Wait a few seconds until the residual vacuum is compensated.
- 6) Now pull the KABEVETTE® together with the Luer cannula out of the vein.
- 7) Disconnect the Luer cannula from the KABEVETTE® only then and dispose it in a sharps disposal container according to the operating instructions. Do not put the protective cannula cap back on!
- 8) Slowly press any air bubbles that may have formed out of the KABEVETTE®, covering the opening with a cellulose swab to absorb any escaping blood.
- 9) Then put the conus cap back on the Luer cone and press it firmly.
- 10) Mix KABEVETTE® carefully and gently through rolling it between the palms 10 times!
- 11) Perform blood gas test within 10 min after collection. If the blood gas test cannot be performed within 10 min after collection, the blood sample must be stored and transported under ice cooling.


Recommendation: Perform blood gas test within 30 min of blood collection if possible.

Storage

Store dry at +5°C to 25°C and protect from direct sunlight and humidity.

All serious incidents related to the product shall be reported to the manufacturer and the responsible authority of the member state where the user and/or patient is registered.

Symbol meaning

	Item number		Sterilization method
	Lot in YYX0000		Protect from sunlight
	Shelf life in YYYY-MM		Protect from rain
	CE-mark + notified body		Do not use if packaging is damaged
	Manufacturer		Follow information documents
	For single use only		Storage temperature
	Do not sterilize again		IVD product
	Single sterile barrier system		

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KABE LABORTECHNIK GmbH