

Bionet tissue embedding cassettes

NEW



Bionet tissue embedding cassettes

Bionet tissue embedding cassettes with dense network for the performance of ultrafine biopsies. Due to the low permeability the insertion of a pad is not necessary; even smallest tissue samples do not penetrate the net. The writing space allows the labelling via a printer or with conventional pens.

Bionet cassettes are available as follows:

Colour	single ¹ with lid	single ² without lid	in a magazine ³ without lid	stacked on a tie ⁴ with lid	stacked on a tie ⁵ without lid
white	054770	054750	054560	054720	054640
light blue	054771	054751	054561	054721	054641
yellow	054772	054752	054562	054722	054642
pink	054773	054753	054563	054723	054643
light green	054774	054754	054564	054724	054644
light orange	054775	054755	054565	054725	054645
brown	054776	054756	054566	054726	054646
red	054777	054757	054567	054727	054647
grey	054778	054758	054568	054728	054648
beige	054779	054759	054569	054729	054649
light purple	054780	054760	054570	054730	054650
turquoise	054781	054761	054571	054731	054651
dark orange	054782	054762	054572	054732	054652
dark purple	054783	054763	054573	054733	054653
pastel yellow	054784	054764	054574	054734	054654
pastel blue	054785	054765	054575	054735	054655
pastel green	054786	054766	054576	054736	054656

Item no.	Description
053831	Lid for tissue embedding cassettes Bionet natural ²

¹ Packing unit: 500 pcs in dispenser box, 3,000 pcs in box

² Packing unit: 500 pcs in bag, 5,000 pcs in box

³ Packing unit: 1,200 pcs in box (16 magazines of 75 pcs),
suitable for printing systems by Thermo Scientific: MicroWriter™ and PrintMate™

⁴ Packing unit: 600 pcs in dispenser box (15 x 40 pcs), 3,600 pcs in box

suitable for printing systems by:

SAKURA: Tissue-Tek AutoWrite

Leica: Automatic printing systems

Vogel: cassette printer

PRIMERA: cassettes printer

⁵ Packing unit: 1,200 pcs in dispenser box (30 x 40 pcs), 7,200 pcs in box

