

FFP2 masks & Medical surgical face masks



The FFP2 masks satisfy the European standard EN 149:2001+A1:2009, are CE-marked and have been tested accordingly by a notified body.

Item No.	Description
*	Respirator FFP2/KN95

Packing unit: on request



The CE-certified medical-surgical face masks manufactured in Germany meet the requirements of ISO 14683 - Medical-Surgical-Face Masks Type II R- and are therefore fully surgery compatible.

Item No.	Description
096001	Medical-surgical face mask type II R

Packing unit: 40 pcs. in dispenser box, 2,760 pcs. in box