

Disposable blades for microtomes



MS200

Universal blade with plasma coating for daily routine: suitable for all types of tissue, especially good for ligament cuts, very good cutting quality even with very thin cuts, suitable for cutting samples embedded in paraffin and plastic.

Extra hardened blade cutting edge, long durability.

Item No.	Description
093981	Microtome disposable blades MS 200



MS300

Very sharp special blade for soft and hard tissue, especially bone tissue, biopsies, mamma, lymph nodes and ligament sections.

Item No.	Description
093982	Microtome disposable blades MS 300



MS33

Blade with a 30° angle allows extra fine cuts on all types of tissue, including skin samples. Short incision time, long life.

Item No.	Description
093983	Microtome disposable blades MS 33



MS24

Blade suitable for hard tissue, but also usable for soft tissue and biopsies, ligament sections and for use on the cryostat.

Item No.	Description
093984	Microtome disposable blades MS 24

Packing unit: each 50 pcs. in box, 10 boxes in dispenser box

