

Test tubes with physiological saline solution & Swabs for throat and Nasopharynx*



The test tubes with physiological saline solution NaCl 0.9% serve as transport vessels for swab samples. The tubes are available in various sizes and designs and are supplied sterile. On request, the stoppers are also available in other colours as well as with personalized (barcode) label.

NaCl-solution Tube (L x Ø) Type of Stopper	NaCl-solution		
	1 ml	1.5 ml	2 ml
75/12 mm, overlapping stopper, white	042500 100 pcs/bag 2,000 pcs/box	-	-
75/13 mm, screw cap, natural	0961 3360 100 pcs/bag 2,000 pcs/box	0961 3363 100 pcs/bag 2,000 pcs/box	-
75/15 mm ¹ , screw cap, natural	0961 3030 100 pcs/bag 1,500 pcs/box	-	0961 3031 100 pcs/bag 1,500 pcs/box
92/15 mm, screw cap, natural	0961 3250 100 pcs/bag 1,000 pcs/box	-	0961 3251 100 pcs/bag 1,000 pcs/box
100/15 mm, overlapping stopper with cams, natural	-	-	042511 100 pcs/bag 1,000 pcs/box
101/16,5 mm, screw cap, natural	0961 1250 100 pcs/bag 1,000 pcs/box	-	0961 1251 100 pcs/bag 1,000 pcs/box

¹ On request, this tube is also available with a raised screw cap. This allows the use of the swabs with predetermined breaking point mentioned below.

The swabs with nylon flock fiber tip are available in two versions for throat and/or nasopharyngeal swab tests. They are individually sterile packed.

Item no.	Description	Packing unit
056507*	Nylon flock fiber swab, nasopharynx, breakpoint 80 mm, tip-Ø 3 mm	50 pcs./bag, 5,000 pcs./box
056509*	Nylon flock fiber swab, throat, breakpoint 80 mm, tip-Ø 5 mm	200 pcs./bag, 5,000 pcs./box