

# Stainless steel sample rack



## Features

The brushed stainless steel sample rack is an ideal companion in workday life of laboratories. It is extremely sturdy, stainless and therefore long-living. Advantages:

- heat-resistant up to 200°C
- autoclavable
- UV-sterilisable
- dishwasher proof
- chemically resistant

## Capacities and dimensions

As standard configuration:

- 2 x 10 bores for vessels with a diameter of Ø 27 mm
- approx. 35.5 cm x 7.5 cm x 3.5 cm (L x W x H)

**Different dimensions and bores available on request!**

## Sample Rack

Item no.	Description
061714	Stainless steel sample rack 2 x 10 bores; Ø 27.5 mm

## Suitable pathology vessels from the KABE LABORTECHNIK portfolio

Item no.	Description	Volume
053650	Pathology vessels with enclosed screw cap	20 ml

Packing unit: 500 pcs in bag, 1,000 pcs in box

Item no.	Description	Volume
042020	Pathology vessels with 7.5 ml formalin 4% buffered	20 ml

Packing unit: 500 pcs in box

