


Urine analysis

URIVETTE®/URIVETTE® BAK: Extraction-, transport- and analysis vessel!



URIVETTE®
(8 ml and 10 ml,
sterile and unsterile)



URIVETTE® BAK
(10 ml, unsterile)



Urine cup with screw
(100 ml, sterile and
unsterile)



Urine collection bottle
(2000 ml and 2500 ml)

URIVETTE® is a disposable urine collection device with scale, screw cap and suction tube.
URIVETTE® BAK offers the stabilisation of the bacterial count in the urine for about 24h.

- Hygienic, clean drawing with suction tube
- Test-strips can be checked in URIVETTE® BAK immediately
- Centrifugal type for the extraction of sediment
(around 5 min with 400-800 g)
- Unbreakable transport vessel
- Available in different sizes, volumes, non-sterile and individual sterile packaging