

Pathology/Histology



Consumables for the treatment of
pathological/histological samples
by KABE LABORTECHNIK

Quality, reliability and consistency since 1977!



Drawing on our long years of experience in medical technology, KABE LABORTECHNIK develops and produces high-quality consumables for daily work in the clinic, medical practice and laboratory. At the same time, our aim is to offer our customers and partners innovative and optimal solutions that make everyday work simpler, safer and more efficient.

Perfection right down to the smallest detail, hygienically immaculate conditions and consistent quality management are the highest priority in manufacturing our products. Although our products are mass-produced, every item is subjected to the same criteria as if it were made to order.

Our product range is as varied as the applications in which our products are put to use. With this brochure, we would like to provide you with a comprehensive overview of our product range in pathology/histology.



Table of contents



Tissue embedding cassettes
(standard, universal, biopsy)

Page 4

Accessories for embedding cassettes

Page 8



Treated test tubes for pathological/
histological analyses

Page 9

Accessories for treated test tubes

Page 9

Pathology vessels with screw cap

Page 11

Colour palette tissue embedding cassettes

	white
	pastel yellow
	yellow
	light orange
	dark orange
	red
	pink
	light purple
	dark purple
	pastel blue
	light blue
	turquoise
	pastel green
	light green
	beige
	brown
	grey



Tissue embedding cassettes Standard



Item no. (without lid)	Colour	Item no. (with lid)
053700	white	053745
053701	light blue	053746
053702	yellow	053747
053703	pink	053748
053704	light green	053749
053705	light orange	053750
053706	brown	053751
053707	red	053752
053708	grey	053753
053709	beige	053754
053710	light purple	053755
053711	turquoise	053756
053712	dark orange	053757
053713	dark purple	053758
053861	pastel yellow	053879
053862	pastel blue	053880
053863	pastel green	053881

Packing unit: without lid – 500 pcs in bag, 5,000 pcs in box
with lid – 500 pcs in dispenser box, 3,000 pcs in box

Tissue embedding cassettes Universal



Item no. (without lid)	Colour	Item no. (with lid)
053715	white	053760
053716	light blue	053761
053717	yellow	053762
053718	pink	053763
053719	light green	053764
053720	light orange	053765
053721	brown	053766
053722	red	053767
053723	grey	053768
053724	beige	053769
053725	light purple	053770
053726	turquoise	053771
053727	dark orange	053772
053728	dark purple	053773
053867	pastel yellow	053885
053868	pastel blue	053886
053869	pastel green	053887

Packing unit: without lid – 500 pcs in bag, 5,000 pcs in box
with lid – 500 pcs in dispenser box, 3,000 pcs in box

Tissue embedding cassettes

- The high-quality material is resistant to solvents that are usually found in the laboratory and guarantees a high level of dimensional stability.
- Labelling is possible on three sides. The front side is slightly sloped and has a special surface texture that ensures long-lasting labelling.
- HB pencils, solvent-resistant markers and commercial labelling machines can be used for labelling.
- Can be supplied without, with separate or with pre-attached hinged lid.

Tissue embedding cassettes Biopsy



Item no. (without lid)	Colour	Item no. (with lid)
053730	white	053775
053731	light blue	053776
053732	yellow	053777
053733	pink	053778
053734	light green	053779
053735	light orange	053780
053736	brown	053781
053737	red	053782
053738	grey	053783
053739	beige	053784
053740	light purple	053785
053741	turquoise	053786
053742	dark orange	053787
053743	dark purple	053788
053873	pastel yellow	053891
053874	pastel blue	053892
053875	pastel green	053893

Packing unit: without lid – 500 pcs in bag, 5,000 pcs in box
with lid – 500 pcs in dispenser box, 3,000 pcs in box

Tissue embedding cassettes stacked




Tissue embedding cassettes stacked in a magazine (without lid)

The embedding cassettes are pre-stacked in a magazine ready for use in cassette printers. One magazine contains 75 embedding cassettes without lids. The embedding cassettes are available in 17 colours in the variants standard, universal and biopsy.

Colour	Item no. Standard	Item no. Universal	Item no. Biopsy
white	054500	054520	054540
light blue	054501	054521	054541
yellow	054502	054522	054542
pink	054503	054523	054543
light green	054504	054524	054544
light orange	054505	054525	054545
brown	054506	054526	054546
red	054507	054527	054547
grey	054508	054528	054548
beige	054509	054529	054549
light purple	054510	054530	054550
turquoise	054511	054531	054551
dark orange	054512	054532	054552
dark purple	054513	054533	054553
pastel yellow	054514	054534	054554
pastel blue	054515	054535	054555
pastel green	054516	054536	054556

Packing unit: 1,200 pcs in box

Compatible lids are available separately, see page 8.



Tissue embedding cassettes stacked on a tie (without and with lid)

The embedding cassettes are pre-stacked on a holding tie ready for use in cassette printers. 40 embedding cassettes without or with lids are threaded onto the tie. The cassettes are available in 17 colours in the variants standard, universal and biopsy.

Without lid

Colour	Item no. Standard	Item no. Universal	Item no. Biopsy
white	054580	054600	054620
light blue	054581	054601	054621
yellow	054582	054602	054622
pink	054583	054603	054623
light green	054584	054604	054624
light orange	054585	054605	054625
brown	054586	054606	054626
red	054587	054607	054627
grey	054588	054608	054628
beige	054589	054609	054629
light purple	054590	054610	054630
turquoise	054591	054611	054631
dark orange	054592	054612	054632
dark purple	054593	054613	054633
pastel yellow	054594	054614	054634
pastel blue	054595	054615	054635
pastel green	054596	054616	054636

Packing unit: 1,200 pcs in dispenser box, 7,200 pcs in box

Compatible lids are available separately, see page 8.

With lid

Colour	Item no. Standard	Item no. Universal	Item no. Biopsy
white	054660	054680	054700
light blue	054661	054681	054701
yellow	054662	054682	054702
pink	054663	054683	054703
light green	054664	054684	054704
light orange	054665	054685	054705
brown	054666	054686	054706
red	054667	054687	054707
grey	054668	054688	054708
beige	054669	054689	054709
light purple	054670	054690	054710
turquoise	054671	054691	054711
dark orange	054672	054692	054712
dark purple	054673	054693	054713
pastel yellow	054674	054694	054714
pastel blue	054675	054695	054715
pastel green	054676	054696	054716

Packing unit: 600 pcs in dispenser box, 3,600 pcs in box

Accessories for embedding cassettes



Plastic lid for embedding cassettes

single-use

Item no.	Description
053800	Lid for embedding cassettes standard
053815	Lid for embedding cassettes universal
053830	Lid for embedding cassettes biopsy

Packing unit: 500 pcs in bag, 5,000 pcs in box



Metal lid for embedding cassettes

re-usable

Item no.	Description
053860	Metal lid for embedding cassettes



Filter paper for embedding cassettes

Filter paper made from fine pored cellulose tissue for holding tissue samples. The filter paper is porous for solvents and fits in all commercial embedding cassettes.

Item no.	Description
094015	Filter paper for embedding cassettes

Packing unit: 50 pads in box (1 pad = 100 sheet)



Foam inserts for embedding cassettes

Foam inserts made from fine pored fabric for holding tissue samples. The foam inserts are porous for solvents and fit in all commercial embedding cassettes.

Item no.	Description
094005	Foam inserts for embedding cassettes blue
094006	Foam inserts for embedding cassettes orange
094016	Foam inserts for embedding cassettes black

Packing unit: 500 pcs in bag

Treated test tubes
for pathological/
histological analyses

2 ml, Ø 11 mm, L = 40 mm



Test tubes 2 ml, prefilled with 4 % formalin solution
as fixative for pathological and histological samples

Phosphate buffered solution pH 7.0 ± 0.3

Item no.	Description
042011	FO 40 L yellow flat bottom 1 ml formalin

Packing unit: 100 pcs in bag, 4,000 pcs in box

The test tube with screwed on stopper is suitable
as a protective vessel, see page 10.

Two test tubes 2 ml fit into the protective vessel.

5 ml, Ø 18 mm, L = 40 mm



Test tubes 5 ml, prefilled with 4 % formalin solution
as fixative for pathological and histological samples

Phosphate buffered solution pH 7.0 ± 0.3

Item no.	Description
042015	FO 18/40 PP L yellow 2 ml formalin

Packing unit: 100 pcs in bag, 1,000 pcs in box

Suitable protective vessel UG 126 S, see page 10.

Two test tubes 5 ml fit into the protective vessel.

Can be supplied with
individual labelling on
request for an additional
charge.



Test tubes prefilled with 4 % formalin solution
as fixative for pathological and histological samples.

Phosphate buffered solution pH 7.0 ± 0.3

Item no.	Description
042020	Pathology vessel 20 ml with 7.5 ml formalin

Packing unit: 500 pcs in box

042022	Pathology vessel 35 ml with 10 ml formalin
--------	--

Packing unit: 500 pcs in box

042023	Pathology vessel 50 ml with 20 ml formalin
--------	--

Packing unit: 100 pcs in box

042021	Pathology vessel 100 ml with 40 ml formalin
--------	---

Packing unit: 100 pcs in box

Accessories for
treated test tubes



Rack for 20 test tubes 2 ml and 5 ml for
pathological and histological analyses

Item no.	Description	Test tube
061701	20/10	2 ml
061711	20/18.5	5 ml



Test tube with screwed on stopper

The test tube with screwed on stopper is suitable as a protective vessel for treated test tube 2 ml, see page 9.

Item no.	Description	Dimensions
0961 3110	N 101/16.5 S	Ø 16.5 mm, L = 101 mm

Packing unit: 1,000 pcs in bag, 1,000 pcs in box



Protective vessels without absorbing inlay

Vessels with transparent screw cap for the transportation of medical and biological specimen.

Item no.	Description	Dimensions
090403	UG 85	Ø 30 mm, L = 85 mm

Packing unit: 500 pcs in bag, 500 pcs in box

090401	UG 114	Ø 44 mm, L = 114 mm
--------	--------	---------------------

Packing unit: 250 pcs in bag, 250 pcs in box

090405	UG 126	Ø 30 mm, L = 126 mm
--------	--------	---------------------

Packing unit: 300 pcs in bag, 300 pcs in box



Protective vessels with absorbing inlay and symbol printing

Vessels with transparent screw cap for the transportation of medical and biological specimen.

- Comply with DIN EN 829
- Protective vessel autoclavable at 121°C/2.1 bar
- Absorbing inlay completely encloses the sample tube
- Optimal protection against external forces

Item no.	Description	Dimensions
090404	UG 85 S	Ø 30 mm, L = 85 mm

Packing unit: 500 pcs in bag, 500 pcs in box

090406	UG 126 S	Ø 30 mm, L = 126 mm
--------	----------	---------------------

Packing unit: 300 pcs in bag, 300 pcs in box



Transport packing for diagnostic specimen

with liquid-tight secondary packing and absorbing inlay.

Comply with ADR-packing instructions P 650 and the packing regulations of Deutsche Post AG.

Item no.	Description
090060	Transport packing for diagnostic specimen

Outer diameter:
230 x 220 x 35 mm
Empty weight: ca. 100 g

Pathology vessels with enclosed screw cap

TWIN-Function:
The pathology vessel
35 ml can be used as
a protective vessel for the
pathology vessel 20 ml.

The pathology vessel
100 ml can be used as
a protective vessel for the
pathology vessel 50 ml.



Item no.	Description
053650	Pathology vessels 20 ml with enclosed screw cap

Packing unit: 500 pcs in bag, 1,000 pcs in box

053651	Pathology vessels 35 ml with enclosed screw cap
053652	Pathology vessels 50 ml with enclosed screw cap

Packing unit: 500 pcs in bag, 500 pcs in box

053653	Pathology vessels 100 ml with enclosed screw cap
--------	--

Packing unit: 250 pcs in bag, 250 pcs in box

KABE LABORTECHNIK – because health is a matter of trust!



In addition to the pathology/histology products,
KABE LABORTECHNIK offers a diverse range of
laboratory materials for the medical laboratory:

- Venous blood collection systems for collecting human medical samples
(Primavette® S and V, KABEVETTE® G and V)
- Venous blood collection systems for collecting veterinary medical samples
(Primavette®, KABEVETTE®)
- Capillary blood collection system for small sample volumes
(Capillary blood collection system GK)
- Plastic and glass capillaries
- Extensive range of accessories for the above blood collection systems
and techniques
- Treated test tubes and collection vessels for blood, urine and
faeces analyses
- Untreated test tubes and reaction vessels, laboratory items made of plastic
- Devices and consumables for blood sugar determination
- Devices and systems for determining the blood sedimentation rate
- Mixers for blood samples and other liquid samples

We would be happy to provide you with our full catalogue.
Or visit us online at
www.kabe-labortechnik.de

KABE LABORTECHNIK GmbH
Jägerhofstraße 17
D-51588 Nümbrecht-Elsenroth
Phone: +49 2293 596
Fax: +49 2293 3823
info@kabe-labortechnik.de
www.kabe-labortechnik.de