Test Tube 15 ml with Formalin


– The perfect solution for secure sample preservation and cost-effective shipping –

With the new test tube 15 ml with formalin, KABE-Labortechnik GmbH is expanding its line of high-quality test tubes for pathology and histology.

Designed to meet the new requirements for letter mail

Starting in January 2027, a **maximum height of 3 cm** and a **maximum weight of 2,000 g** will apply to **basic products for national letters**; this affects the “Large Letter” and “Maxi Letter” formats. From January 2027, **mailings that are larger and/or heavier** than the new “Maxi Letter” format must be sent using **products from the parcel portfolio**.*

The new sample container was specifically designed for these changed shipping requirements. With its compact design, it enables space-saving and therefore cost-effective shipping while offering a reliable solution for daily use in pathology and histology.



Item no.	Description	Closure	Packaging
042850	Test tube 15 ml with 4 % formaldehyde solution	screw lid	50 pieces in dispenser box

Item no.	Description
061772	Rack for Test tubes 15 ml, 2 rows á 10 bores made of stainless steel

Your benefits at a glance:

- ✓ Pre-dosed phosphate-buffered 4% formaldehyde solution, pH 7.0 ± 0.3 as a ready-to-use fixative
- ✓ Guaranteed leak-tightness in accordance with DIN EN ISO 6717:2021 and packaging instruction P650
- ✓ Yellow screw-on lid for internal identification of 4 % formaldehyde solution at KABE
- ✓ Standard label included
- ✓ Optional: Custom labeling according to customer specifications
- ✓ Additional lid colors available on request
- ✓ Proven KABE quality

* Mail & Parcel Germany | Product Changes for Letters | March 2026

KABE-Labortechnik GmbH • Jägerhofstraße 17 • D-51588 Nümbrecht-Elsenroth
Phone: + 49 2293 9132-0 • Fax: + 49 2293 9132-7030 • E-Mail: info@kabe-labortechnik.de • www.kabe-labortechnik.de